ABDO BELT S™
ABDOBELTLINE.COM

MANUAL

APOLYNE® - ABDO BELT S™

ENGLISH

Technical Specifications

- Power: 12 W
- Voltage: 100-240 V
- Frequency: 50-60 Hz

Included

- 1 Abdo Belt S™ belt
- 1 extension for the belt
- 1 cable with plug for a 240 V socket
- 1 cable with plug adaptor for 12 V socket
- 1 150 cm measuring tape

Instructions

1. Connect the belt to the mains
2. Place the Abdo Belt S™ on the area of the body that you wish to treat, securing it comfortably without overtightening it. You can use the Abdo Belt S™ on your shoulders, abdomen, back, waist, hips, thighs and calves.
3. To access the automatic vibration function, press the

power button. The AUTO/MANUAL indicator light will turn red. It will start at the lowest power level, gradually rising and falling in intensity.

4. To access the manual vibration mode, turn the device on and press the SET button, the AUTO/MANUAL indicator light will turn green. It will start at the lowest intensity. Press the + and - buttons to increase or decrease the vibration intensity which will be displayed on the SPEED LEVEL indicator.
5. In the sauna function you can choose the time duration (15, 30 or 50 minutes) and the temperature (35 °C, 40 °C or 45 °C). To access this mode, press the TIMER button and the lowest duration and temperature values will automatically be selected. Press the TIMER button several times to change the duration, which will appear on the display. To select the temperature, press the TEMPE button until the display shows the desired temperature.
6. The vibration and sauna functions can be used simultaneously or individually, in both automatic and manual modes.

Warnings

Turn the Abdo Belt ST™ off immediately if you notice any discomfort. Adjust the speed and temperature of the device based on your capabilities. If you suffer from any health problems, consult a doctor before using the device. The Abdo Belt ST™ is not recommended for use by pregnant women. Only use this device on the body

parts indicated in the instruction manual. If it is used for a long period of time, the overheating protection may become active and this will turn off the appliance immediately; when the temperature decreases you can reuse it again normally. Unplug the device when it is not in use. Do not use the Abdo Belt STM in high-temperature environments, near any source of water or in high levels of humidity. Clean it with a slightly damp cloth. Do not submerge this device in liquid. Do not let the product or its parts come into contact with hot surfaces. Keep the Abdo Belt STM in a dry and dust-free place. If the appliance, its cable or its components seem damaged or they are not functioning properly, do not use them or do not attempt to repair them; only a qualified technician should repair them. Keep it out of the reach of children and pets. This product is not household waste. Please follow the environmental regulations in place for its recycling.

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