[WWW.ABDOKING.COM](http://WWW.ABDOKING.COM)

**ABDO KING**

---

# **ENGLISH**

## **Operation**

1. Install 2 x AAA batteries (included) in the battery compartment.
2. Attach the controller to the front part of the belt.
3. Adjust the belt with the hook-and-loop fasteners. The front and back parts must be centred.
4. Remove the belt and attach the electrodes to it using the press studs.
5. Apply plenty of gel to the electrodes.
6. To turn on the unit, press the on/off key.
7. Use the key with three rectangles to set the different options: mode (arrow), intensity (bars) and timer (clock).
8. Select the desired mode:
  1. Long contractions
  2. Very fast pulses
  3. Tapping with rest periods
  4. Mixed exercise with rest periods
  5. Exercise with increasing intensity
  6. Mixed exercise
  7. Front and back tapping
  8. Full exercise
9. Press the key with three rectangles and then the Front or Back keys to select the intensity (10 levels).
10. Select the timer and set it by pressing the Front or Back keys. The default time is 10 min and it can be adjusted up to 45 min.
11. To flip the display, press the side button.
12. To turn on or off the sound, press the display key with an image of a speaker and keep it pressed for 2 seconds.
13. To lock or unlock the keys, press the display key with an image of a key and keep it pressed for 3 seconds.
14. To switch off the unit, press the on/off key.

## **Cautions**

1. The adjustable belt has different sizes: S (61-73.5 cm), M

- (81-91.5 cm) and L (99-112 cm). If you need a bigger size, extend the belt using the extension band with hook-and-loop fasteners.
2. For the unit to work, the electrodes must have gel and be in contact with the skin. Use water-based gels only, never oil-based gels. If you try to turn on the unit without wearing it or without gel, an alarm will sound.
  3. Do not bend the electrodes. When they are worn out, replace them with new electrodes (not included).
  4. Do not place electrodes on the navel.
  5. The unit can be used on the waist, chest, legs and arms, but never on the neck or head.
  6. During the first week of training, only use modes 1 to 4.
  7. The unit will switch off automatically after 10 min when not in use.
  8. Switch off the unit before removing it.
  9. Adapt the training intensity to your physical condition.
  10. Do not wash the belt or use fabric softener on it. To clean the belt, use a slightly damp cloth. Do not immerse the belt in liquids.
  11. The unit is not suitable for people with a pacemaker or heart problems, epilepsy, multiple sclerosis, phlebitis, inflamed tissues, varicose veins in its later stage or cancer. Do not use it on period days or if you have given birth recently.
  12. The unit is not suitable for use by children unless they are supervised by a responsible adult.
  13. The belt is not a medical device.
  14. Store the unit when not in use.

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